

Amon

Single-Use Micro-Grasping
Forceps Needle for
Intrascleral Fixation Techniques



Multifunctional Precision



My forceps needle is a significant innovation in ophthalmic surgery that enhances the efficiency and precision of complex procedures, particularly in IOL fixation.

It was designed to provide intraocular access like a needle while also functioning as a forceps to grasp tissue, suture material or haptics in a single maneuver.




Michael Amon

Prof. Michael Amon (Austria)



Original size

The Amon 27G Single-Use Micro-Grasping Forceps Needle



Challenge

Intrascleral haptic fixation of IOLs is a complex procedure that requires high precision. The use of conventional needles involves a secondary instrument inserting the haptic into its lumen, with difficulty increasing for the trailing haptic. Disengaging is also a potential risk when using conventional needles.

Solution

A combination of needle-shaped probe with internal micro-grasping forceps. The instrument provides intraocular access, grasps and secures haptics into its lumen and externalizes them, minimizing the risk of disengagement. All in one instrument, in one continuous manoeuvre.

The Amon Forceps Needle Offers a Broad Range of Applications

- Scleral fixation of multiple IOL designs
- Iris fixation of IOLs
- Fixation of capsular tension rings
- Fixation of capsular segments
- Gore-Tex and Polypropylene sutures




Watch the following video for an in-vitro and in-vivo demonstration by Prof. Michael Amon (Austria)
<https://bit.ly/41j9ZFb>

Features of the Amon Forceps Needle


Beveled needle-shaped tip with excellent piercing properties for sclerotomies and scleral tunnels

Micro-grasping forceps facilitate the handling of IOL haptics and suture materials




Secure grasping ability prevents unexpected disengaging of IOL haptics or suture materials


Create
scleral tunnel




Rotate for
optimal position



Grasp
IOL haptic



Hold and
externalize safely



In conclusion, the instrument's 2-in-1 design increases efficiency by enabling one-step maneuvers. This saves time in the operating room and minimizes the risk of complications.



The Yamane Flanged Intrasccleral Haptic Fixation (FIHF) Technique with Amon Forceps Needle


The Amon Forceps Needle was developed to address certain limitations encountered with flanged intrasccleral haptic fixation (FIHF) techniques such as described by Dr. Shin Yamane (Japan).

In comparison to conventional needles, the Amon Forceps Needle prevents slippage of IOL haptics during externalization. As the forceps needle enhances control and precision, it further simplifies and standardizes the Yamane technique with its numerous advantages, if used with the Yamane Double Needle Stabilizer by Geuder:

- Simplified localization of the sclerotomy positions
- Control of puncture angles to form intrasccleral tunnels
- Standardized insertion angles
- Fixation of the eye during needle insertion



Yamane Double Needle Stabilizer, original size (details on page 11)




The Yamane Double Needle Stabilizer provides precise localization of the entry points as well as insertion and inclination angles for the needle when creating the scleral tunnels, in a reproducible way. The addition of the Amon Forceps Needle and its innovating and simplifying features, are a step forward towards standardization of intrasccleral fixation techniques.

Step by Step

Yamane Technique with Amon Forceps Needle


1. Stabilization

The Yamane Double Needle Stabilizer features globe fixation as an improvement over the single fixation point that tissue forceps can achieve.




2. Precise insertion

The platforms offer visual confirmation of the entry points at the respective distance from limbus and the tangential angle of insertion. As well as the inclination angle in respect to the iris plane.




3. Simplified manoeuvre

The Amon Forceps Needle is suitable to work in combination with the stabilizer, blending the precision and simplifying features of both surgical tools.



"The Forceps Needle is a Brilliant Innovation."

What Ophthalmologists Say



"The Forceps Needle
is a brilliant innovation."

Robert H. Osher, MD

The Video Journal of Cataract,
Refractive and Glaucoma Surgery



Watch the
video here



"It made my life easier by making
the surgery more efficient and easier
for me and that translates into
easier recovery for the patient."

Brandon D. Ayres, MD

www.youtube.com/@brandonayres



Watch the
video here



"... provides a simplification of
Dr. Yamane's procedure making
it easier for surgeons to perform
this technique for ISHF."

Sadeer B. Hannush, MD

2024 American Academy of Ophthalmology



Watch the
video here

Prof. Michael Amon in Focus

**41ST ANNUAL
ASCRS**
film festival

Runner-up price for Prof. Amon in the category new techniques:
Utilizing needle function and grasping function
with a new instrument for different scleral fixation
techniques of IOLs: the Forceps Needle


[41st Annual ASCRS Film Festival](#)

Category of new techniques

 **ESCRS**
2024

At the European Society of Cataract and Refractive
Surgeons meeting, Prof. Michael Amon's **"Saloon Door"**
technique wins **ESCRS Video Award 2024**. The technique
is used during an open-sky IOL exchange with the help
of the Forceps Needle.

[ESCRS video awards 2024](#)



Recommended literature:

Flanging Techniques in Anterior Segment
Surgery – Mastering Aphakia,
IOL-exchange and IOL-refixation




[Prof. Michael Amon](#)

<https://bit.ly/4h68aBj>

Specifications and Order Information

G-34502

AMON
Single-Use Forceps Needle
for intrascleral IOL fixation
27 gauge / 0.4 mm
6 pcs. per box, sterile
for use with JAPAN MODEL Handle G-38240



G-38240

JAPAN MODEL
Handle
for capsule scissors and forceps
interchangeable tips
6 mm, titanium
overall length 90 mm



Recommended Accessories

G-19828

Single-Use Low Temperature Battery Cautery
fine tip (1.5 mm), tip length 28 mm,
max. temp. 600 °C
10 pcs. per box, sterile



Recommended Accessories for Yamane Technique

The Amon Forceps Needle can be used in combination with the Yamane Double Needle Stabilizer to overcome limitations of conventional cannulas or needles used for externalizing IOL haptics.

G-31497


YAMANE

Double Needle Stabilizer

Fixation ring for intrascleral IOL fixation
with double needle technique
for 27 gauge needles
overall length 110 mm



The Snyder-Osher IOL explantation scissors and forceps facilitate the removal of IOLs through small incisions avoiding induced astigmatism.



SNYDER-OSHER

IOL Scissors

For cutting acrylic and silicone foldable IOLs
lower notched blade to create IOL traction
18 gauge / 1.2 mm shaft
overall length 150 mm
short model, 3.5 mm blades
long model, 4.5 mm blades



G-32998

G-32996

G-32997

SNYDER-OSHER

IOL Forceps

to hold the IOL while cutting it inside
the capsular bag
2 mm long jaws with 0.5 mm serrated tips
22 gauge / 0.7 mm lightly curved shaft
overall length 99 mm



