


## STERILIZATION CONTAINERS AND ADAPTERS



MICROSTOP <sup>®</sup> .....	<b>249</b>
STEEL CASES .....	<b>252</b>
TRAYS & RACKS FOR STEEL CASES .....	<b>256</b>
SILICONE STRIPS & MATS .....	<b>257</b>
SMALL CASES .....	<b>258</b>
STERISAFE <sup>®</sup> DURO .....	<b>259</b>
ADAPTERS .....	<b>264</b>



**G-39745**

**MICROSTOP<sup>®</sup>  
CONTAINER LOW**

for one tray  
outside dimensions: 310 x 189 x 70 mm  
inside dimensions: 283 x 177 x 50 mm



**G-39750**

**MICROSTOP<sup>®</sup>  
CONTAINER HIGH**

for two trays  
outside dimensions: 310 x 189 x 90 mm  
inside dimensions: 283 x 177 x 70 mm



**G-39783**

**MICROSTOP<sup>®</sup>  
TRAY MINI**

29 mm, low  
size: 277 x 172 x 29 mm

without lid

**G-39786**

with lid



**G-39785**

**MICROSTOP<sup>®</sup>  
TRAY MINI**

38 mm, high  
size: 277 x 172 x 38 mm

without lid


**G-39784**

with lid




MICROSTOP<sup>®</sup>

**G-39787**      **MICROSTOP®**  
**TRAY LID MINI**  
 size: 277 x 172 mm  
 for G-39783 and G-39785




**G-39782**      **MICROSTOP®**  
**RETAINER ELEMENT**  
 for MICROSTOP® tray mini  
 nap Style, 120 mm  
 separable into 6 single elements  
**G-39789**      29 mm, low  
 38 mm, high



**G-39788**      **MICROSTOP®**  
**CLIP FOR RETAINER ELEMENTS**  
 G-39782 and G-39789  
 10 pcs. per bag



**G-39791**      **MICROSTOP®**  
**ID LABEL**  
 for tray, aluminium  
 customized labeling  
 silver-colored



**G-39741**      **MICROSTOP®**  
**LABEL**  
 for lid, customized labeling



**G-39746**      **MICROSTOP®**  
**ID LABEL**  
 for front, customized labeling




**MICROSTOP<sup>®</sup>  
LOGISTIC FRAME**  
for lid to use with labels

- G-39753** color: red
- G-39754** color: blue
- G-39755** color: green




- G-39762** **MICROSTOP<sup>®</sup>  
SEAL**  
with indicator  
1000 pcs. per box



- G-39764** **MICROSTOP<sup>®</sup>  
TAG**  
with indicator  
1000 pcs. per box



- G-39680** **PROTOCOL CARDS**  
for ERGO-SAFE tray  
320 pcs. per box



**G-19822**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 220 x 165 x 55 mm
- textile filter in lid and base
- silicone seal
- perforated insert without edge
- perforated silicone tray G-32785
- 2 labels with clamp



**G-19823**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 220 x 165 x 75 mm
- textile filter in lid and base
- silicone seal
- perforated insert without edge
- perforated silicone tray G-32785
- 2 labels with clamp



**G-28151**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
- silicone seal
- perforated insert without edge
- 3 silicone burl strips
- 2 silicone paddings in lid
- 2 labels with clamp



**G-28150**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert with edge, subdivided, with perforated silicone plate and perforated silicone tray
- silicone insert, small, with perforated silicone tray
- 2 labels with clamp



**G-19876****STERILIZATION CASE**

- for steam sterilization  
outside dimensions: 300 x 220 x 55 mm
- textile filter in lid and base
  - silicone seal
  - perforated insert without edge
  - 3 silicone burl strips
  - 2 silicone paddings in lid
  - 2 labels with clamp

**G-19856****STERILIZATION CASE**

- for steam sterilization  
outside dimensions: 300 x 220 x 55 mm
- paper filter with quick-clasp in lid and base
  - silicone seal in the lid
  - without insert and profile strips
  - 2 labels with clamp

**G-19843****STERILIZATION CASE**

- for steam sterilization  
outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
  - silicone seal
  - perforated insert with edge
  - perforated silicone plate
  - perforated insert without edge
  - 3 silicone burl strips
  - 2 silicone paddings in lid
  - 2 labels with clamp



**G-19842**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert with edge and perforated silicone plate
- perforated insert with edge perforated silicone plate and silicone instrument tray
- 2 labels with clamp



**G-19877**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 400 x 290 x 55 mm
- textile filter in lid and base
- silicone seal
- perforated insert without edge
- 4 silicone burl strips
- 3 silicone paddings in lid
- 2 labels with clamp



**G-19892**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 400 x 290 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- 2 perforated inserts without edge with 4 silicone nap strips each
- 3 silicone profile strips at the base of the upper perforated insert
- 3 silicone profile strips in the lid
- 2 labels with clamp






**G-19894**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 400 x 290 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert without edge with 4 silicone nap strips
- perforated insert with edge and silicone nap plate
- 3 silicone profile strips in the lid
- 2 labels with clamp



**G-19890**

**STERILIZATION CASE**

- for steam sterilization
- outside dimensions: 495 x 270 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert with edge 4 silicone nap strips and lateral compartment (240 x 110 x 40 mm),
- 3 silicone profile strips in the lid
- 2 labels with clamp



**PERMANENT TEXTILE FILTER**

for steam sterilization cases

- G-31650** for G-19822, G-19823
- G-31651** for G-19841, G-19842, G-19843, G-19876, G-28150, G-28151
- G-31652** 6 perforations, for G-19877, G-19892, G-19894
- G-31656** for G-19877
- G-31661** for G-19890



**G-31649**

**INDICATOR**


for steam sterilization, slim  
1000 pcs. per box




**G-32675**     **SILICONE TRAY**  
with bottom  
outside dimensions: 160 x 85 x 20 mm




**G-32660**     **STERILIZATION RACK**  
for 12 cannulas, silicone  
outside dimensions: 85 x 56 x 8 mm



**G-19870**     **STERILIZATION RACK**  
for 10 instruments with handle, silicone  
outside dimensions: 160 x 85 x 20 mm



**G-32785**     **SILICONE PERFORATED TRAY**  
outside dimensions: 194 x 135 x 24 mm




**G-39650**      **STERILIZATION CASE**  
for steam sterilization, plastic  
for 6 vitreoretinal instruments  
outside dimensions: 210 x 130 x 26 mm



**G-32676**      **SILICONE NAP STRIP BLUE**  
340 x 30 mm




**G-32677**      **SILICONE STRIP BLUE**  
350 mm long




**SILICONE BURL MAT**  
perforated

**G-19866**      240 x 230 mm  
**G-19865**      300 x 150 mm  
**G-19867**      470 x 245 mm



**G-19899**      **SILICONE PLATE**  
perforated, 270 x 185 mm  
for hot air sterilization cases  
G-19841, G-19842, G-19843,  
G-19846, G-19847



SILICONE STRIPS & SILICONE MATS


**G-31657 SAFETY SEAL**  
for steam sterilization cases  
100 pcs. per box



**G-31659 SEAL**  
for hot-air sterilization cases  
100 pcs. per box



**G-32695 SILICONE CLOSURE**  
17 cm long



**G-33080 STERILIZATION CASE**  
for steam sterilization, plastic, for 5 instruments with handle  
outside dimensions: 160 x 85 x 20 mm



**G-33090 STERILIZATION CASE**  
for one instrument or diamond knife, plastic  
outside dimensions: 160 x 60 x 20 mm



**G-33100 STERILIZATION CASE**  
for one vitreoretinal instrument, plastic  
outside dimensions: 195 x 60 x 25 mm



**G-36090**      **STERILIZATION CASE**  
plastic, for 6 vitreoretinal lenses



**G-39527**      **STERISAFE<sup>®</sup> DURO A4 PLUS**  
container for steam, gas (EO, FO) and plasma  
sterilization, without trays and strips  
outside dimensions: 290 x 190 x 70 mm  
inside dimensions: 265 x 165 x 60 mm



**G-39526**      **STERISAFE<sup>®</sup> DURO A4 PLUS COMPLETE SET**  
container for steam, gas (EO, FO) and plasma  
sterilization, including:  
1 x Sterisafe<sup>®</sup> wire basket A4-30  
5 x Sterisafe<sup>®</sup> silicone retainer strips 230 mm  
1 x Sterisafe<sup>®</sup> metal tray A4-30  
outside dimensions: 290 x 190 x 70 mm  
inside dimensions: 265 x 165 x 60 mm

**G-39525**      **STERISAFE<sup>®</sup> DURO A4**  
container for steam, gas (EO, FO) and plasma  
sterilization  
outside dimensions: 290 x 190 x 60 mm  
inside dimensions: 265 x 165 x 49 mm



**G-39528**      **STERISAFE<sup>®</sup> DURO A6**  
container for steam, gas (EO, FO) and plasma  
sterilization  
outside dimensions: 190 x 110 x 60 mm  
inside dimensions: 165 x 92 x 49 mm



STERISAFE<sup>®</sup> DURO


**G-39543 STERISAFE® DURO WIRE BASKET A4-30**

outside dimensions: 260 x 163 x 30 mm  
inside dimensions: 249 x 157 x 23 mm




**G-39541 STERISAFE® DURO METAL TRAY A4-15**

outside dimensions: 260 x 163 x 15 mm  
inside dimensions: 249 x 157 x 11 mm



**G-39542 STERISAFE® DURO METAL TRAY A4-30**

outside dimensions: 260 x 163 x 30 mm  
inside dimensions: 249 x 157 x 26 mm




**G-39549 STERISAFE® DURO WIRE BASKET A6-30**

outside dimensions: 162 x 90 x 30 mm  
inside dimensions: 156 x 80 x 23 mm



**G-39548 STERISAFE® DURO METAL TRAY A6-15**

outside dimensions: 162 x 90 x 15 mm  
inside dimensions: 156 x 80 x 11 mm



**G-39545**      **STERISAFE<sup>®</sup> SILICONE RETAINER STRIP 75**  
for Sterisafe<sup>®</sup> DURO A4 and A6  
length: 75 mm



**G-39546**      **STERISAFE<sup>®</sup> SILICONE RETAINER STRIP 145**  
for Sterisafe<sup>®</sup> DURO A4 and A6  
length: 145 mm



**G-39547**      **STERISAFE<sup>®</sup> SILICONE RETAINER STRIP 230**  
for Sterisafe<sup>®</sup> DURO A4 and A6  
length: 230 mm




**G-39501**      **STERISAFE<sup>®</sup> FILTER**  
A4 10 pcs. per box  
**G-39511**      A6 10 pcs. per box



**G-39530**      **STERISAFE<sup>®</sup> FILTER TAG EXPIRATION DATE**  
10 pcs. per box



**G-39535**      **SEAL / LABEL**  
seal with color indicator and blank label  
2000 pcs. per box


Steri. - Datum/Date	
Name	
	Dampf (braun)  Steam (brown)

STERISAFE<sup>®</sup> DURO


**G-39540**      **SECURITY SEAL**  
100 pcs. per box




**G-39502**      **STERISAFE® SILICONE BURL MAT**  
A4  
**G-39512**      A6



**G-39503**      **STERISAFE® SILICONE GRID INSERT**  
A4  
**G-39513**      A6




**G-39504**      **STERISAFE® METAL GRID INSERT**  
A4  
**G-39514**      A6





**G-39505**     **STERISAFE<sup>®</sup> METAL GRID INSERT**  
with 2 silicone burl strips



**G-31670**     **MELAG METAL GRID INSERT**  
with 2 silicone nap strips  
outside dimensions: 370 x 147 x 50 mm  
for MELAG sterilizing case

