

**BAUSCH + LOMB**

# RETINA

Product Catalog



# Table of Contents



**Stellaris Elite**  
Vision Enhancement System

VITRECTOMY **4**

INSTRUMENTATION **22**

LASER PROBES **36**

ILLUMINATION **40**

OTHER PRODUCTS **46**



# ENHANCED RESEARCH AND DEVELOPMENT

Because retinal surgery is always changing, we are dedicated to staying ahead of the curve. With a unique combination of talent and resources, we constructed our integrated research and development team to do just that. By preparing for what will come, rather than reacting to what has already happened, we are striving to change the way you approach surgery.





vit  
VITRECTOMY



## Hardware + Accessories

### Stellaris Elite™

BL15455

- Combined System with integrated 532nm laser, includes 2 lamp modules (Xenon and/or Mercury/Xenon).

BL14455

- Combined System without integrated 532nm laser, includes 2 lamp modules (Xenon and/or Mercury/Xenon).

### Zero Level Bottle Hanger

BL4363 Dual Bottle Hanger



### Integrated Foot Pedal

BL2295



### Auxiliary footpedal for LIO use

BL2296

### Integrated Foot Pedal Accessories

- BL4390 Foot Control Battery
- BL4391 Foot Control Wall Charger
- BL4393 Foot Control Charging Cradle
- BL4394 Foot Control Power Cable



## Fragmentation Handpiece

BL3270



## Stellaris® Phaco Handpiece

BL3170

- Needle and sleeve not included




## Mayo Tray Adapter

BL4364




## Laser filters

- BL3236 Zeiss Fixed ESF
- BL3240 Zeiss 2 position filter ESF
- BL3245 Moeller Wedel Fixed ESF
- BL3241 Leica 2 position filter ESF
- BL3247 Leica 4-way filter
- BL3237 Leica Fixed ESF
- BL3239 Laser filter ESF Mount Topcon
- BL3242 Adapter for fitting of 2 position filters



## Xenon Lamp Replacement

BL3320



## Fragmentation Needle 15°

20G BL5275

- Single use



## Fragmentation Needle 15°


20G MVS1600A

- Reusable



## Mercury Lamp Replacement

BL3321



## Fragmentation Needle 30°

23G BL5277

- Single use



## Laser Safety Goggles

BL3235

## Interlock bypass connector kit

BL3244

## Laser Upgrade kit

BL2272

## LIO

BL2298

## LIO Replacement Bulb

BL3246

## Laser On/Off Keys (1 set)

BL3234


# Stellaris Elite™ Posterior Surgical Packs

## 20G Posterior Pack

SE5420M | Mid-Field  
SE5420W | Widefield

### Pack Includes:

- A 20G Vitrectomy Cutter
- B 20G Light pipe
- C 20G MVR knife
- D Scleral plug package - 19 and 20




## 23G Posterior Pack

SE5423MV | Mid-Field  
SE5423WW | Widefield

### Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



## 25G Posterior Pack

SE5425MV | Mid-Field  
SE5425WW | Widefield

### Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



### All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Pre-piercing spike assembly
- 9 Stop cock



1



2



3



4



5



6



7



8



9




# Stellaris Elite™ Combined Surgical Packs

## 20G Combined Pack

SE5520M | Mid-Field  
SE5520W | Widefield

### Pack Includes:

- A** 20G Vitrectomy Cutter
- B** 20G Light pipe
- C** 20G MVR knife
- D** Scleral plug package - 19 and 20




## 23G Combined Pack

SE5523MV | Mid-Field  
SE5523WV | Widefield

### Pack Includes:

- A** 23G Vitrectomy Cutter
- B** 23G Light pipe
- C** 23G Valved ESA System
- D** 23G Scleral plugs



## 25G Combined Pack

SE5525MV | Mid-Field  
SE5525WV | Widefield

### Pack Includes:

- A** 25G Vitrectomy Cutter
- B** 25G Light pipe
- C** 25G Valved ESA System
- D** 25G Scleral plugs



### All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Irrigation sleeves x 2
- 9 Needle wrench
- 10 Pre-piercing spike assembly
- 11 Stop cock
- 12 Test chamber




## Stellaris®PC Posterior X Surgical Packs with Bi-Blade®

### 27G Posterior Pack

BL5427MVXB

**Pack Includes:**

- A** 27G Bi-Blade® Vitrectomy Cutter
- B** 27G Mid-field light pipe
- C** 27G Valved ESA System




## Stellaris Elite™ Posterior Surgical Packs with Bi-Blade®

### 25G Posterior Pack

25G SE5425MVB | Midfield  
25G SE5425WVB | Widefield

**Pack Includes:**

- A** 25G Bi-Blade® Vitrectomy Cutter
- B** 25G Light pipe
- C** 25G Valved ESA System
- D** 25G Scleral plugs




### 27G Posterior Pack

27G SE5427MVB

**Pack Includes:**

- A** 27G Bi-Blade® Vitrectomy Cutter
- B** 27G Mid-field light pipe
- C** 27G Valved ESA System



**All the above packs come with the items below.**

- 1** Cassette with AFI tube assembly
- 2** Second Aspiration Line
- 3** Irrigation Line
- 4** Drapes for Mayo tray, screen and remote control
- 5** Cassette stopper
- 6** Priming tray
- 7** Male luer coupler
- 8** Pre-piercing spike assembly
- 9** Stop cock



1



2



3



4



5



6



7



8



9


## Stellaris®PC Combined X Surgical Packs with Bi-Blade®

### 27G Combined Pack with Bi-Blade®

BL5527MVXB

#### Pack Includes:

- 27G Bi-Blade® Vitrectomy Cutter
- 27G Mid-field light pipe




## Stellaris Elite™ Combined Surgical Packs with Bi-Blade®

### 25G Combined Pack

25G SE5525MVB | Midfield  
25G SE5525WVB | Widefield

#### Pack Includes:

- A** 25G Bi-Blade® Vitrectomy Cutter
- B** 25G Light pipe
- C** 25G Valved ESA System
- D** 25G Scleral plugs




### 27G Combined Pack

27G SE5527MVB

#### Pack Includes:

- A** 27G Bi-Blade® Vitrectomy Cutter
- B** 27G Mid-field light pipe
- C** 27G Valved ESA System



### All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Irrigation sleeves x 2
- 9 Needle wrench
- 10 Pre-piercing spike assembly
- 11 Stop cock
- 12 Test chamber




# Stellaris®PC Posterior X Surgical Packs with AFI (Air Forced Infusion)

## 20G Posterior Pack

BL5420WX

### Pack Includes:

- A 20G Vitrectomy Cutter
- B 20G Wide angle light pipe
- C 20G MVR knife
- D Scleral plug package - 19 and 20




## 23G Posterior Pack

BL5423WVX

### Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



## 25G Posterior Pack

BL5425WVX

### Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



### All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Pre-piercing spike assembly
- 9 Stop cock



1



2



3




4



5



6



7



8



9



# Stellaris®PC Posterior X Surgical Packs (Gravity)

## 20G Posterior Pack

BL5220WX

### Pack Includes:

- A 20G Vitrectomy Cutter
- B 20G Wide angle light pipe
- C 20G MVR knife
- D Scleral plug package - 19 and 20





## 23G Posterior Pack

BL5223WVX

### Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs





## 25G Posterior Pack

BL5225WVX


### Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



### All the above packs come with the items below.

- 1 Cassette with gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Stop cock




# Stellaris®PC Combined X Surgical Packs with AFI (Air Forced Infusion)

## 20G Combined Pack

BL5520WX

### Pack Includes:

- A** 20G Vitrectomy Cutter
- B** 20G Wide angle light pipe
- C** 20G MVR knife
- D** Scleral plug package - 19 and 20




## 23G Combined Pack

BL5523WVX

### Pack Includes:

- A** 23G Vitrectomy Cutter
- B** 23G Wide angle light pipe
- C** 23G Valved ESA System
- D** 23G Scleral plugs



## 25G Combined Pack

BL5525WVX

### Pack Includes:

- A** 25G Vitrectomy Cutter
- B** 25G Wide angle light pipe
- C** 25G Valved ESA System
- D** 25G Scleral plugs



### All the above packs come with the items below.

- 1** Cassette with AFI tube assembly
- 2** Second Aspiration Line
- 3** Irrigation Line
- 4** Drapes for Mayo tray, screen and remote control
- 5** Cassette stopper
- 6** Priming tray
- 7** Male luer coupler
- 8** Irrigation sleeves x 2
- 9** Needle wrench
- 10** Pre-piercing spike assembly
- 11** Stop cock
- 12** Test chamber





# Stellaris®PC Combined X Surgical Packs (Gravity)

## 20G Combined Pack

BL5320WX

### Pack Includes:

- A 20G Vitrectomy Cutter
- B 20G Wide angle light pipe
- C 20G MVR knife
- D Scleral plug package - 19 and 20





## 23G Combined Pack

BL5323WVX

### Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs





## 25G Combined Pack

BL5325WVX


### Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



### All the above packs come with the items below.

- 1 Cassette with combined gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming tray
- 7 Male luer coupler
- 8 Irrigation sleeves x 2
- 9 Needle wrench
- 10 Stop cock
- 11 Test chamber



## Stellaris®PC

BL14304

- Combined System with integrated 532nm laser
- Includes 2 lamp modules: Xenon and/or Mercury/Xenon

BL14334

- Combined System without integrated 532nm laser
- Includes 2 lamp modules: Xenon and/or Mercury/Xenon




## Valved Trocars

### 23G ESA Kit-Valved

BL5281


- Includes infusion cannula, tubing and plugs
- (Kit of 3)



### 25G ESA Kit-Valved

BL5282


- Includes infusion cannula, tubing and plugs
- (Kit of 3)



### 27G ESA Kit-Valved

BL5240

- Includes infusion cannula, and tubing
- (Kit of 3)



### Single Trocar-Valved

23G BL5283

25G BL5285

27G BL5241

### Replacement Valves

23G BL5284

25G BL5286

ESA valved cannula system designed to deliver remarkable performance and exceptional surgeon control



## Vitrectomy Cutters for Stellaris®PC

20G VIT Cutter  
BL5620



23G VIT Cutter  
BL5623S



25G VIT Cutter  
BL5625S



## Vitrectomy Cutters for Stellaris Elite™

20G VIT Cutter  
BL5626 - ELITE



23G VIT Cutter  
SE5623 - ELITE



25G VIT Cutter  
SE5625 - ELITE



## Vitrectomy Handle Extension

BL3394  
• Box of 12



## Bi-Blade® Vitrectomy Cutters

### For Stellaris Elite™

25G SE5625B  
27G SE5627B

### For Stellaris®PC

27G BL5627S


**bi-Blade®**



## MVR Knives

### LaserEdge® Plus

19G LP7519 A: 23mm B: 1.4mm  
20G LP7520 A: 23mm B: 1.2mm  
23G LP7523 A: 23mm B: 0.8mm  
25G LP7525 A: 23mm B: 0.6mm



## Vitreoretinal Tamponades

### Universal Gas Kit

TR906061




### High Flow Viscous Fluid Injection Cannulas

20G BL7620  
23G BL7623  
25G BL7625



### Viscous Fluid Injection Adapter

BL7610



### Cylinder Stand

TR9066



### Adato® Sil-Oil 5000 centiStokes

ES5000S

- 10 mL
- Syringe to be used with CX5720



### Silicone Oil Collection Chamber

30.90

- Silicone tip allows for removal of 1000-5000 cSt silicone oil through an existing instrument cannula
- Extrusion controlled via foot pedal
- 15mL collection chamber volume
- Sterile, 6/box



### Standard Regulator

TR9085



### Viscous Fluid Extraction Kits

CX5710



### 10 Dose Cylinder of SF<sub>6</sub>

TR9088



### Viscous Fluid Injection Kits for Oxane Range

CX5720



### 10 Dose Cylinder of C<sub>3</sub>F<sub>8</sub>

TR9089



## VersaVIT™ 2.0 Standard Procedure Packs

### Standard Pack w/VersaVIT™ Midfield

20G 70020S  
23G 70023S  
25G 70025S



#### Pack Includes:

- Vitrectomy Cutter
- Cassette
- 3cc Syringe
- 3-Way Stopcock
- Screen Drape
- Adapter Kit
- Light Pipe
- Valved Port Entry System (except 20G)
- Soft Tip Cannula (except 20G)
- Adjustable Gas-Pressured Infusion Tubing
- Extrusion Tubing
- MVR Blade (20G only)

### Standard Pack w/PHOTON Widefield

20G 70020SP  
23G 70023SP  
25G 70025SP




#### Pack Includes:

- Vitrectomy Cutter
- Cassette
- 3cc Syringe
- 3-Way Stopcock
- Screen Drape
- Adapter Kit
- Light Pipe
- Valved Port Entry System (except 20G)
- Soft Tip Cannula (except 20G)
- Adjustable Gas-Pressured Infusion Tubing
- Extrusion Tubing
- MVR Blade (20G only)

## VersaVIT™ 2.0 Core Essential Vitrectomy + Core Essential Biopsy Packs

### Core Essentials™ Pack

23G 70023E  
25G 70025E




#### Pack Includes:

- Vitrectomy Cutter
- Cassette
- 3cc Syringe
- 3-Way Stopcock
- Screen Drape
- Adapter Kit

### Core Essentials™ Biopsy Pack

23G 70123B  
25G 70125B



#### Pack Includes:

- Vitrectomy Cutter
- Cassette
- 3cc Syringe
- 3-Way Stopcock
- Screen Drape
- Adapter Kit
- Biopsy Vial Kit

## VersaVIT™ 2.0 Vitrectomy Accessories

### Owner's Manual

VWITOME

### Adapter Hose — Y Connector

VWITAHOSE-Y

### Adapter Hose Extended

VWITAHOSE-E

### Replacement Fuse

VWITFUSE-1.0 1 Amp  
VWITFUSE-2.0 2 Amp




## VersaVIT™ Vitrectomy Accessories

### Adjustable Gas-Pressurized Infusion Tubing Set

17565-10

- Includes 1 prespike and 2 stopcocks
- 6/box



### Viscous Fluid Injection System

VFIHP/VFIMOD


- Linear control of injection pressure through a reusable, titanium, hand piece
- Maximum injection pressure of 100 psi (as dictated by facility's medical gas supply regulator)
- Custom sterilization tray included
- For use with Synergetics® VFI packs



### IV Standard Gravity Administration Set

17570

- 6/box



### Viscous Fluid Injection Pack

19020 Polyimide VFI cannulas

19030 Stainless Steel VFI Cannulas

- Provides the surgeon the consumables necessary to perform viscous fluid injection: syringe, plunger rod, piston, syringe tip cap, injection airline and VFI cannulas (20ga, 23ga and 25ga)
- For use with Synergetics® VFI System
- Sterile, 6/box



### Fluid-Air Exchange Tubing

17680

- 6/box



### Vitrectomy Cutter

23G 70933

25G 70935

- For use with VersaVIT 2.0™
- 6,000 CPM

### Extrusion Tubing Set

17580

- 6/box



---

### Silicone Oil Collection Chamber

30.90

- Silicone tip allows for removal of 1000-5000 cSt silicone oil through an existing instrument cannula
- Extrusion controlled via foot pedal
- 15mL collection chamber volume
- Sterile, 6/box



---

### Port Entry System w/ Valves

23G 58.23STHV

25G 58.25STHV

- Set of 3



---

### VersaVIT™ Midfield Endoilluminator

20G 56.65

23G 56.65.23

25G 56.65.25

27G 56.65.27



---

### Dual Infusion Kit

BL5280



---

### Port Entry System

20G 58.20-6

- Polyimide valved cannula

- Single packed



---

### DigiFlow™ Tubing

D4600A





A close-up, black and white photograph of a microfluidic instrument. The image shows several white, cylindrical components, likely syringes or reservoirs, arranged in a row. A thin, clear tube is connected to one of the components. The background is a soft, out-of-focus white. A teal square frame is overlaid on the center of the image, containing the text 'in' and 'INSTRUMENTATION'.

in  
INSTRUMENTATION

## Handle Designs

### Pinnacle 360°™

PIN

Pinnacle 360° is our line of fully disposable, hand-held retinal surgical instruments. As its name suggests, the Pinnacle 360° product line provides high quality and value. With Pinnacle 360°, surgeons can expect consistent performance in a variety of applications.

- Offers wide variety of disposable tip options
- Fully disposable
- 360° actuation to approach the retina from any angle
- Ergonomic design to reduce fatigue
- Color coded handles for gauge identification
- 20ga, 23ga, 25ga, 27ga, Foreign Body



### Syntrifugal™ Handle and Disposable Tips

D00.01SL, D00.01SS

The Syntrifugal line provides surgeons with a reusable titanium handle with a disposable tip. Surgeons enjoy consistent tactile feel with the freedom to choose from a comprehensive selection of disposable tips.

- Large and small handle options
- Handle design that allows the tip to rotate 360°
- 20ga, 23ga, 25ga, 27ga



### Premium

Our fully reusable Premium line delivers consistency and performance. Surgeons enjoy the familiar comfort and reliability of knowing their tools will provide exceptional performance day in and day out.

- Large and small handle options
- Add "S" to the end of catalog number to order small handle
- Blackened stainless steel tips to reduce glare
- Etched serial numbers for quick and easy identification
- Linear closure
- Microserrations™ offered for select tip styles
- 18ga, 19ga, 20ga, 23ga, 25ga



PIN Pinnacle 360°™  
D Disposable Tip used with Syntrifugal Handle  
# ONLY Premium reusable instrument



# ILM

---

## Eckardt Forceps

20G 39.08PIN; D39.08; 39.28  
23G 39.08.23PIN; D39.08.23; 39.28.23  
25G 39.08.25PIN; D39.08.25; 39.28.25  
27G 39.08.27PIN



---

## End Grasping Forceps

20G 09.12PIN  
23G 09.12.23PIN  
25G 09.12.25PIN



---

## Fine Tip Eckardt Forceps

20G 09.10  
23G 09.10.23PIN; D39.10.23; 09.10.23  
25G 09.10.25PIN; D39.10.25; 09.10.25



---

## Lambert ILM Forceps

20G 09.26



---

## Microserrated™ End Grasping Forceps

23G 39.12.23PIN  
25G 39.12.25PIN



---

## Morris/Witherspoon Fullview ILM Forceps

20G D09.40; 09.50  
25G D09.40.25; 09.50.25



---

## Pincer End Grasping Forceps

20G 11.28  
23G 11.28.23



---

## Pointed Micro Forceps

20G 11.23  
23G 11.23.23



---

## Tano Asymmetrical Micro Forceps

20G 09.04PIN; D09.04; 39.23  
23G 39.03.23PIN; D39.03.23; 39.23.23  
25G 39.03.25PIN; D39.03.25; 39.23.25  
27G 39.03.27PIN; D39.03.27



## ILM

### Microserrated™ Snub Nosed Forceps

23G 39.18.23PIN



## Dual Purpose

### Apte Corset Forceps

23G 39.07.23PIN, 39.27.23  
25G 39.07.25PIN, 39.27.25



### Tano Asymmetrical Peeling Forceps

20G 39.24



### Microserrated™ Tewari Omnicep Forceps

23G 39.15.23PIN, D39.15.23  
25G 39.15.25PIN, D39.15.25  
27G 39.15.27PIN



## Heavy Membrane

### Pyramid Tip Forceps

20G 11.04PIN, 11.24  
23G 11.04.23PIN  
25G 11.04.25PIN



## Foreign Body

### BB Removal Forceps

18G 11.40PIN, 11.40



### Serrated Jaw Forceps

20G 11.06PIN, 11.26  
23G 11.06.23PIN,  
25G 11.06.25PIN



### Foreign Body Forceps

19G 11.20PIN, 11.20



### IO Foreign Body Magnet

18G 12.21PIN, 12.21



PIN Pinnacle 360™

D Disposable Tip used with Syntrifugal Handle

# ONLY Premium reusable instrument

# Scissors

---

## Horizontal Scissors Curved

20G 23.76PIN, D23.76  
23G 23.76.23PIN, D23.76.23  
25G 23.76.25PIN, D23.76.25



---

## Morris Vertical Scissors

20G 10.24PIN, 10.24  
23G 10.24.23PIN, 10.24.23  
25G 10.24.25PIN, 10.24.25



---

## Horizontal Scissors Medium Curved

20G 23.26  
23G 23.26.23  
25G 23.26.25



---

## Straight Scissors

20G 23.29  
23G 23.29.23



---

## Horizontal Scissors Small Curved

20G 23.27



---

## Vertical Scissors

20G 10.21



---

## Microserrated™ Horizontal Scissors Medium Curved

20G 23.56  
23G 23.56.23  
25G 23.56.25



## Pick Forceps

### Ducournau Serrated Forceps

20G D29.02, 29.22  
23G 29.22.23



## Subretinal

### Thomas Horizontal Subretinal Forceps

20G 39.21



## Myopic

### Fine Tip Eckardt Forceps

23G 09.10.23PIN-M  
25G 09.10.25PIN-M  
• 37mm



## Pick Forceps

20G 09.05PIN, 09.25  
23G 09.05.23PIN, 09.25.23  
25G 09.05.25PIN, 09.25.25



### Thomas Horizontal Subretinal Scissors 135°

20G 23.23



### Microserrated™ Tewari Omnicap Forceps

23G 39.15.23PIN-M  
25G 39.15.25PIN-M  
• 37mm



## Vertical Forceps

23G 11.41.23



### Thomas Vertical Subretinal Forceps

20G 21.21



## Illuminated Disposables

### Irrigating Illuminated Capone Pick

20G 56.44P



### Illuminated Capone/Rizzo Pick

23G 56.45.23P



### Illuminated Membrane Rake

23G 56.46.23P



## Scleral Markers

### Awh Sutureless

51.20

- Marked at 3mm and 4mm



### Awh

51.35


- Marked at 3mm and 4mm



### Awh Horseshoe

51.45

- Marked at 3mm and 4mm
- Specifically designed to stabilize cannula



### Chow

51.25

- Marked at 3mm and 4mm



## Reusable Picks + Spatulas

---

### Rice ILM Elevator 127°

- 20G 20.02
- 23G 20.02.23
- 25G 20.02.25
- 42ga
- .10mm at tip



---

### Shiraga Pick, 110°

- 20G 20.31
- 33ga
- .20mm at tip



---

### Awh Serrated Pick 135°

- 20G 20.41
- 23G 20.41.23
- 25G 20.41.25
- 30ga
- .30mm at tip



---

### Membrane Peeler, Spherical Tip 110°

- 20G 28.04
- 30ga
- .30mm at tip



---

### Michels' Hook 90°

- 20G 28.21
- 36ga
- .20mm at tip



## Reusable Picks + Spatulas

---

### Thomas Subretinal Perforator/Spatula 36ga 131°

20G 20.01

- 36ga
- .10mm at distal tip




---

### Lewis 33ga Subretinal Perforator 135°

20G 27.01

- 36ga
- .20mm at tip



---

### Irrigating Spatula 120°

20G 28.30


- 24ga
- .55mm at tip



---

### Rounded Tip Spatula 140°

20G 28.07



---

### Thomas Blunt 36ga Subretinal Perforator Spatula 131°

20G 20.30

- 36ga
- .20mm at tip




---

### Membrane Rake, Flat Tip 140°

20G 28.05 Flat tip

20G 28.06 Curved tip



---

### Lewis 33ga Subretinal Perforator 90°

20G 26.01

- 36ga
- .20mm at tip



---

### Stellaris®PC Single Use VitreoRetinal Membrane Pick

23G SUK01

25G SUK02



# Reusable Cannulas + Needles

## Reflux Handle

31.00

- For use with products 30.01–30.09 and 30.76



## Universal Soft Tip

30.10

30.10.23

- Reusable Soft Tip with extendable silicone tip



## Blunt Extrusion Needle

20G 30.01



## Soft Tip Cannula

20G 30.08



## Replacement Reflux Reservoir

31.01 2/box,  
31.02 10/box



## Tapered Blunt Extrusion Needle

20G 30.02

- Tapers to 27ga
- .40mm tip



## Resnick Subretinal Fluid Drainage Needle

25G 30.09

- 2.5mm long 1mm bevel



## Thomas Subretinal, 30ga Needle

20G 02.03

- 30ga
- .30mm tip



## Brush Cannula Tip, Straight

20G 30.06



## Metal Locking Luer-Hub Needle

- With products in this section



PIN Pinnacle 360™™

D Disposable Tip used with Syntrifugal Handle

# ONLY Premium reusable instrument



---

### Chang Perfluorocarbon Cannula

- 20G 08.01
- 25ga
  - .50mm tip



---

### Thomas Subretinal, 33ga Cannula

- 20G 13.01
- 33ga
  - .20mm tip



---

### Lewis Subretinal, 33ga TPA Cannula 90°

- 20G 15.01
- 33ga
  - .20mm tip



---

### Thomas Subretinal, 33ga Sharp Bent Cannula

- 20G 12.01
- 33ga
  - .20mm tip



---

### Lewis Subretinal, 33ga TPA Cannula 135°

- 20G 14.01
- 33ga
  - .20mm tip



---

## Disposable Backflush Instruments

---

### Disposable Backflush Handle



---

### Blunt Tip

- 20G D31.01
- 5/box



---

### Brush Tip

- 20G D31.06
- 5/box



---

### Soft Tip

- 20G D31.08  
23G D31.08.23  
25G D31.08.25
- 5/box



## Disposable Cannulas + Needles

### Rigid Micro Injection Cannula, Straight

20G 12.02, 12.03

- .10mm at tip
- Black, disposable handle
- 5/box



### Rigid Micro Injection Cannula, Angled

20G 12.12, 12.13

- .10mm at tip
- Black, disposable handle
- 5/box



### De Juan/Awh Subretinal Injection

25G 12.03.25

- 41ga tip
- .10mm at tip
- 5/box



### Soft Tip Cannula

23G 30.08.23

25G 30.08.25

27G 30.08.27

- 6/box



### Slotted Soft Tip Cannula

23G 30.18.23

25G 30.18.25

- 6/box



### Vitreoretinal Soft Tip

20G 30.76


- For use with 31.00 handle
- 31mm tubing length
- 3mm flexible silicone
- 10/box



### Transconjunctival Chow Silicone Oil Injection Cannula

25G 30.80.25

- 6/box



### Viscous Fluid Injection Pack

19020 Polyimide VFI cannulas

19030 Stainless Steel VFI Cannulas

- Provides the surgeon the consumables necessary to perform viscous fluid injection: syringe, plunger rod, piston, syringe tip cap, injection airline and VFI cannulas (20ga, 23ga and 25ga)
- For use with Synergetics® VFI System
- Sterile, 6/box



## Instrument Trays

### Small

28.25

- 6" x 2" x .75"
- 15.2cm x 5.1cm x 1.9cm
- Holds 2 picks or 1 instrument



### Medium

23102SYN


- 10" x 6" x .75"
- 25.4cm x 15.2cm x 1.9cm
- Holds 10 picks or 6 instruments



### Large

23101SYN


- 15" x 10" x .75"
- 38.1cm x 25.4cm x 1.9cm
- Holds 25 instruments



### Diamond Knife

23103SYN

- 7.5" x 2.6" x .75"
- 19.1cm x 3.8cm x 1.9cm
- Holds 1 diamond knife



## Tip Protectors

---

### Tip Protector, Handled Picks

34.71

- Fits Synergetics® handled picks



---

### Tip Protector for Picks, Rakes and Spatulas

34.72

- Fits Synergetics® handled picks



---

### Tip Protector, Split, Universal

34.74

- Fits Synergetics® Large and Small Squeeze handles



---

### Apex™ Tip Protector

34.78

- Doubles as cleaning luer and tip protector for the APEX™ handle



## Luers + Plugs

---

### Cleaning Luer

34.80

- For Synergetics® premium instruments



---

### Disposable Scleral Plug

19G 12.41

20G 12.40





A close-up, black and white photograph of a white, cylindrical laser probe. The probe is positioned horizontally, with a dark cable extending from its right side. The background is a light, neutral color, possibly a surface or a wall. A teal-colored square frame is overlaid on the center of the image, containing the text 'las' in a large, lowercase, sans-serif font, and 'LASER PROBES' in a smaller, uppercase, sans-serif font below it.

las  
LASER PROBES

## Directional™ Laser Probes

---

### Directional™ Laser Probe

20G 55.26  
23G 55.26.23  
25G 55.26.25, 55.26.25S  
27G 55.26.27



## Straight Laser Probes

---

### Straight Laser Probe

20G 55.21  
23G 55.21.23  
25G 55.21.25S  
27G 55.21.27



## Curved Laser Probes

---

### Flexible Tapered Laser Probe

23G 55.25.23  
25G 55.25.25



### Extendable Directional™ Laser Probe

20G 55.26E  
23G 55.26.23E  
25G 55.26.25E



### TruCurve™ Curved Laser Probe, 37° Curve

20G 55.23  
23G 55.23.23  
25G 55.23.25  
27G 55.23.27



### Inverted Directional™ Laser Probe

20G 55.36E  
23G 55.36.23E  
25G 55.36.25E



### Directional™ II Laser Probe

20G 55.27  
23G 55.27.23  
25G 55.27.25  
27G 55.27.27



DISPOSABLES STERILE, 12/BOX

## Illuminated Laser Probes

### Curved Fixed Extended Illuminated Laser Probe, 37° Curve

20G 55.69P  
23G 55.69.23P  
25G 55.69.25P



### Fixed Extended Illuminated Laser Probe

23G 55.62.23P  
25G 55.62.25P



### Flexible Tapered Illuminated Laser Probe

23G 55.70.23P  
25G 55.70.25P



### Illuminated Directional™ Laser Probe

20G 55.47P  
23G 55.47.23P  
25G 55.47.25P



### Illuminated Directional™ II Laser Probe

23G 55.48.23P  
25G 55.48.25P



## Stellaris® Laser Disposables

### Standard Laser Probe

20G BL5290 STD. Laser PROBE Straight  
23G BL5293 STD. Laser PROBE Straight  
25G BL5295 STD. Laser PROBE Straight



### Curved Standard Laser Probe

20G BL5290C STD. Laser PROBE Curved  
23G BL5293C STD. Laser PROBE Curved  
25G BL5295C STD. Laser PROBE Curved



### Illuminated Laser Probe

20G BL5290L ILL. Laser PROBE Straight  
23G BL5293L ILL. Laser PROBE Straight  
25G BL5295L ILL. Laser PROBE Straight



### Curved Illuminated Laser Probe

20G BL5290LC ILL. Laser PROBE Curved  
23G BL5293LC ILL. Laser PROBE Curved  
25G BL5295LC ILL. Laser PROBE Curved



### Standard Aspirating Laser Probe

20G BL5290ASPH Hard Tip asp. Laser PROBE  
23G BL5293ASPH Hard Tip asp. Laser PROBE



### Soft-tip Aspirating Laser Probe

20G BL5290ASP Soft Tip asp. Laser PROBE S  
23G BL5293ASP Soft Tip asp. Laser PROBE S



DISPOSABLES STERILE, 12/BOX







ILLUMINATION

## Light Source Accessories

---

### Replacement Bulb

PLAMP PHOTON  
P2-PLAMP PHOTON II

## Illumination Adapters

---

### Adapter

57.27 for use with DORC Brightstar  
57.30 for use with Alcon Constellation  
57.40 for use with Bausch + Lomb Stellaris®PC and Stellaris Elite™



## Endoilluminators

---

### Widefield Diffusion

20G 56.21P  
23G 56.21.23P  
25G 56.21.25PS



### Corona Shielded Widefield

20G 56.12P  
23G 56.12.23P  
25G 56.12.25P



### Midfield

20G 56.22P  
23G 56.22.23P  
25G 56.22.25P



### Vivid

20G 56.20P  
23G 56.20.23P  
25G 56.20.25P



### Focal

20G 56.02P  
23G 56.02.23P  
25G 56.02.25PS  
27G 56.02.27P



## Chandeliers

---

### Awh Vivid


25G 56.54.25P

- Trocar/cannula included



### Awh

25G 56.50.25P



### Oshima Vivid

27G 56.54.27P

- Trocar/cannula included



### Oshima Dual

29G 56.50.29P

- Trocar/cannula included



## Illuminated Infusion Chandeliers

---

### Standard

20G 56.30P  
23G 56.30.23P  
25G 56.30.25P



### High-Flow

25G 56.34.25P



## Multifunctional Illumination

---

### Chow Membrane Pick

23G 56.07.23P  
25G 56.07.25PS



### High-Flow Irrigating Endoilluminator

19G 56.14



### Illuminated Bipolar Cautery

25G 56.00.25P



### Aspirating Endoilluminator

23G 56.24.23P  
25G 56.24.25P



## Stellaris® Illumination Disposables

### Focal

20G BL5720  
23G BL5723  
25G BL5725



### Awh Vivid Chandelier

25G BL5235S



### Midfield Endoilluminator

23G BL5823M  
25G BL5825M  
27G BL5827M



### Oshima Vivid Chandelier

27G BL5237S



### Diffuse Widefield Endoilluminator

20G BL5820D  
23G BL5823D  
25G BL5825D



### Adjustable Chandelier

23G BL5234  
25G BL5235

- Trocar/cannula ordered separately




### Wide Angle Endoilluminator

20G BL5820  
23G BL5823  
25G BL5825







op

OTHER PRODUCTS

## DDMS™ Diamond Dusted Membrane Scrapers

---

### Tano

20G 20.04  
23G 20.04.23  
25G 20.04.25S  
27G 20.04.27

- 5/box



### Stiff Tano

20G 20.07

- 5/box



### Retractable Tano

20G 20.16  
23G 20.16.23  
25G 20.16.25

- 3/box—20ga
- 5/box—23ga and 25ga



### Adjustable Tano with Pick

20G 20.18

- 3/box



## Disposable Infusion Cannulas

---

### 6mm Tip

23G 30.32.23

25G 30.32.25

- For severe vitreous hemorrhaging
- For use with Synergetics® steel cannula sets (sold separately)
- Sterile, 6/box



### Infusion Cannula and Infusion Tubing

20G MVS1038A 4mm

23G MVS1069 6mm reusable

25G MVS1070 6mm reusable





## Reusable Infusion Cannulas

---

### 2.5mm Tip

20G 30.20

- Ideal aphakic length
- 1 each



### 4mm Tip

20G 30.21

- All purpose length
- 1 each



### 6mm Tip

20G 30.22

- Ideal aphakic length
- 1 each



## Bipolar

---

### Bipolar Catery

23G 53.00.23

25G 53.00.25

27G 53.00.27

- Sterile, 6/box



## Knives

---

### Opremcak BRVO

20G 52.00

- Super-thin knife, for sheathotomy
- Sterile, 5/box



### Opremcak CRVO

20G 52.20

- Designed for relaxing scleral ring, cribriform plate, and sclera of optic disk during radial optic neurotomy
- Sterile, 5/box



## Iris Retractors

### Reusable

40.01-2 2/box

40.01-6 6/box

- 5 retractors/set in case one is lost or damaged
- Autoclavable up to 3 uses



### Disposable

40.02-2 Sterile, 2/box

40.02-6 Sterile, 6/box

- 5 retractors/set in case one is lost or damaged



## Visco Dissector

### Hubbard Visco Dissector

25G 12.06.25 28ga tip

27G 12.06.27 31ga tip

- Injects viscoelastic for delaminating large membranes from the retina
- Sterile, 6/box



## Stellaris® Bipolar Instruments + Accessories

### Bipolar Eraser

- 23G E7930 Tapered, Lemos type connector
- 23G E7931 Tapered, 2-pin connector
- 25G E7924 Tapered, 2-pin connector
- 25G E7925 Blunt, 2-pin connector



### Bipolar Cord with Banana Plugs

- S2050B
- Reusable



### Banana Plug Adapter

- CX9404
- Reusable



### Bipolar Cord with Banana Plugs

- S2070C
- Single use



### Bipolar Cord

- CX9400
- Reusable



### Bipolar Cautery

- 23G 53.00.23
- 25G 53.00.25
- 27G 53.00.27
- Fully disposable intraocular cautery.
- Single piece design includes handpiece and cord.
- Used with CX9404
- Sterile, 6/box



### Bipolar Cord, Lemo Connector

- CX9430
- Reusable



## I-Pack™ Components

### I-PACK™

Everything you need for office injections in one complete, sterilized I-PACK. The I-PACK will reduce ordering time, reduce and organize your inventory, and reduce preparation time since everything you need is conveniently packaged together. The handy plastic tray also serves as a reservoir that you can use for your povidone-iodine solution. Simply pour the solution into the compartment with the speculum, where it can be used for both your speculum and for dipping cotton swabs. The two additional compartments can be used for discarded vials, swabs, and needles, making cleanup quick and easy.



### I-PACK™ [Sterile, 10/box]

#### 18060 COMPONENTS

- 1mL Tuberculin Syringe
- 1 Pair Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

#### 18061 COMPONENTS

- 2 1mL Tuberculin Syringes
- 1 Pair Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")

#### 18064 COMPONENTS

- 2 1mL Tuberculin Syringes
- 1 Pair Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")
- 1 Disposable Forceps, 4"
- 4 4"x 4" Gauze Pads

#### 18067 COMPONENTS

- 1 30ga 1/2" Needle
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

#### 18069 COMPONENTS

- 1 Pair Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

#### 18069W COMPONENTS

- 1 Pair Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 6 Cotton Swabs 6"
- 1 Wire Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

## Lens Rings

### Awh Sutureless Lens Ring

41.00

- Used to stabilize a vitrectomy lens
- 3 sizes/pack
- Sterile, 10/box



## Disposable Lenses [Sterile, 10/box]

### Flat

VFD

- Field of view: 36°
- Image magnification: 1.0x



### 30° Prism

V30PD

- Field of view: 33° (offset 30°)
- Image magnification: 1.0x



### Bi-Concave

VBCD

- Field of view: 25°
- Image magnification: 0.8x



### Flat, SSV

VFLATSSVD

- Field of view: 30°
- Image magnification: .92x



### Magnifying

VMD

- Field of view: 30°
- Image magnification: 1.5x



### Widefield

VWFD

- Field of view: 48°
- Image magnification: 0.5x



## Pediatric Vitrectomy

### Infusing Instrument Cannulas

58.2023-3MM

- 3mm shaft length
- Stainless steel cannulas
- Allows for two port PPV via 20ga sclerotomy and 23ga instrumentation



### Port Entry System

58.23-3MM

- 3mm shaft length
- Available in 23ga
- Stainless steel cannulas



### Valved Port Entry System

58.23-3MMV

- 3mm shaft length
- Available in 23ga
- Stainless steel cannulas with removable valve cap



## Pediatric Instrumentation

### Vertical Forceps

23G 11.41.23

- Dual function use as a forceps and spatula
- Available in reusable Premium handle configuration



### ILM Microserrated™ Snub Nosed Forceps

23G 39.18.23PIN

- Blunt tip profile, combined with tip cleavage, allows tissue to be engaged in a fashion analogous to two rolling pins moving toward each other
- Available in disposable Pinnacle 360™ handle configuration



### Irrigating Spatula, 120°

20G 28.30

- 24ga (.55mm) at tip
- .55mm at tip
- Reusable



### Irrigating Trese Spatula

23G 28.30.23

- Hand-crafted tip geometry is ideal for dissection
- Glare-free finish
- Standard luer lock connection
- Reusable
- Non-cannula compatible design allows for a longer platform



### Irrigating Trese Spatula

23G D28.30.23

- Cannula compatible design allows for dissection without the need for scleral sutures
- Glare free finish
- Standard luer lock connection
- Sterile, single use, 6/box



### Vertical Scissors

20G 10.24PIN

23G 10.24.23PIN

25G 10.24.25PIN

- Curved heel minimizes inadvertent contact with retina and conforms to eye
- Available in disposable Pinnacle 360™ handle configuration



## Pediatric Illumination

---

### Capone/Rizzo Pick

23G 56.45.23P

- Midfield light pattern
- Semi-malleable pick may be adjusted intraocularly
- Glare-free finish
- Available in all PHOTON configurations



### 2-Port Vitrectomy High-Flow Irrigating Endoilluminator

20G 56.14

- Widefield illumination
- Available in non-PHOTON configuration
- Sterile, 12/box



### Illuminated Membrane Rake

23G 56.46.23P

- Midfield light pattern
- Identify and lift edge of epiretinal membranes
- Glare-free finish
- Available in all PHOTON configurations



### Irrigating Illuminated Capone Pick

20G 56.44P

- Midfield light pattern
- Semi-malleable pick may be adjusted intraocularly
- Glare-free finish
- Available in all PHOTON configurations



# Index

## By Part Number

02.03.....	<b>32</b>	15.01.....	<b>33</b>	30.18.25.....	<b>34</b>	55.26.23E.....	<b>38</b>
08.01.....	<b>33</b>	17565-10.....	<b>20</b>	30.20.....	<b>49</b>	55.26.25.....	<b>38</b>
09.04PIN.....	<b>25</b>	17570.....	<b>20</b>	30.21.....	<b>49</b>	55.26.25E.....	<b>38</b>
09.05PIN.....	<b>28</b>	17580.....	<b>20</b>	30.22.....	<b>49</b>	55.26.25S.....	<b>38</b>
09.05.23PIN.....	<b>28</b>	17680.....	<b>20</b>	30.32.23.....	<b>48</b>	55.26.27.....	<b>38</b>
09.05.25PIN.....	<b>28</b>	18060.....	<b>52</b>	30.32.25.....	<b>48</b>	55.27.....	<b>38</b>
09.10.....	<b>25</b>	18061.....	<b>52</b>	30.76.....	<b>34</b>	55.27.23.....	<b>38</b>
09.10.23.....	<b>25</b>	18064.....	<b>52</b>	30.80.25.....	<b>34</b>	55.27.25.....	<b>38</b>
09.10.23PIN.....	<b>25</b>	18067.....	<b>52</b>	30.90.....	<b>18, 21</b>	55.27.27.....	<b>38</b>
09.10.23PIN-M.....	<b>28</b>	18069.....	<b>52</b>	31.00.....	<b>32</b>	55.36E.....	<b>38</b>
09.10.25.....	<b>25</b>	18069W.....	<b>52</b>	31.01.....	<b>32</b>	55.36.23E.....	<b>38</b>
09.10.25PIN.....	<b>25</b>	19020.....	<b>20, 34</b>	31.02.....	<b>32</b>	55.36.25E.....	<b>38</b>
09.10.25PIN-M.....	<b>28</b>	19030.....	<b>20, 34</b>	34.71.....	<b>35</b>	55.47P.....	<b>39</b>
09.12PIN.....	<b>25</b>	20.01.....	<b>31</b>	34.72.....	<b>35</b>	55.47.23P.....	<b>39</b>
09.12.23PIN.....	<b>25</b>	20.02.....	<b>30</b>	34.74.....	<b>35</b>	55.47.25P.....	<b>39</b>
09.12.25PIN.....	<b>25</b>	20.02.23.....	<b>30</b>	34.78.....	<b>35</b>	55.48.23P.....	<b>39</b>
09.25.....	<b>28</b>	20.02.25.....	<b>30</b>	34.80.....	<b>35</b>	55.48.25P.....	<b>39</b>
09.25.23.....	<b>28</b>	20.04.....	<b>48</b>	39.03.23PIN.....	<b>25</b>	55.62.23P.....	<b>39</b>
09.25.25.....	<b>28</b>	20.04.23.....	<b>48</b>	39.03.25PIN.....	<b>25</b>	55.62.25P.....	<b>39</b>
09.26.....	<b>25</b>	20.04.25S.....	<b>48</b>	39.03.27PIN.....	<b>25</b>	55.69P.....	<b>39</b>
09.50.....	<b>25</b>	20.04.27.....	<b>48</b>	39.07.23PIN.....	<b>26</b>	55.69.23P.....	<b>39</b>
09.50.25.....	<b>25</b>	20.07.....	<b>48</b>	39.07.25PIN.....	<b>26</b>	55.69.25P.....	<b>39</b>
10.21.....	<b>27</b>	20.16.....	<b>48</b>	39.08PIN.....	<b>25</b>	55.70.23P.....	<b>39</b>
10.24.....	<b>27</b>	20.16.23.....	<b>48</b>	39.08.23PIN.....	<b>25</b>	55.70.25P.....	<b>39</b>
10.24PIN.....	<b>27, 54</b>	20.16.25.....	<b>48</b>	39.08.25PIN.....	<b>25</b>	56.00.25P.....	<b>44</b>
10.24.23.....	<b>27</b>	20.18.....	<b>48</b>	39.08.27PIN.....	<b>25</b>	56.02P.....	<b>42</b>
10.24.23PIN.....	<b>27, 54</b>	20.30.....	<b>31</b>	39.12.23PIN.....	<b>25</b>	56.02.23P.....	<b>42</b>
10.24.25.....	<b>27</b>	20.31.....	<b>30</b>	39.12.25PIN.....	<b>25</b>	56.02.25PS.....	<b>42</b>
10.24.25PIN.....	<b>27, 54</b>	20.41.....	<b>30</b>	39.15.23PIN.....	<b>26</b>	56.02.27P.....	<b>42</b>
11.04PIN.....	<b>26</b>	20.41.23.....	<b>30</b>	39.15.23PIN-M.....	<b>28</b>	56.07.23P.....	<b>44</b>
11.04.23PIN.....	<b>26</b>	20.41.25.....	<b>30</b>	39.15.25PIN.....	<b>26</b>	56.07.25PS.....	<b>44</b>
11.04.25PIN.....	<b>26</b>	21.21.....	<b>28</b>	39.15.25PIN-M.....	<b>28</b>	56.12P.....	<b>42</b>
11.06PIN.....	<b>26</b>	23.23.....	<b>28</b>	39.15.27PIN.....	<b>26</b>	56.12.23P.....	<b>42</b>
11.06.23PIN.....	<b>26</b>	23.26.....	<b>27</b>	39.18.23PIN.....	<b>26, 54</b>	56.12.25P.....	<b>42</b>
11.06.25PIN.....	<b>26</b>	23.26.23.....	<b>27</b>	39.21.....	<b>28</b>	56.14.....	<b>44, 55</b>
11.20.....	<b>26</b>	23.26.25.....	<b>27</b>	39.23.....	<b>25</b>	56.20P.....	<b>42</b>
11.20PIN.....	<b>26</b>	23.27.....	<b>27</b>	39.23.23.....	<b>25</b>	56.20.23P.....	<b>42</b>
11.23.....	<b>25</b>	23.29.....	<b>27</b>	39.23.25.....	<b>25</b>	56.20.25P.....	<b>42</b>
11.23.23.....	<b>25</b>	23.29.23.....	<b>27</b>	39.24.....	<b>26</b>	56.21P.....	<b>42</b>
11.24.....	<b>26</b>	23.56.....	<b>27</b>	39.27.23.....	<b>26</b>	56.21.23P.....	<b>42</b>
11.26.....	<b>26</b>	23.56.23.....	<b>27</b>	39.27.25.....	<b>26</b>	56.21.25PS.....	<b>42</b>
11.28.....	<b>25</b>	23.56.25.....	<b>27</b>	39.28.....	<b>25</b>	56.22P.....	<b>42</b>
11.28.23.....	<b>25</b>	23.76PIN.....	<b>27</b>	39.28.23.....	<b>25</b>	56.22.23P.....	<b>42</b>
11.40.....	<b>26</b>	23.76.23PIN.....	<b>27</b>	39.28.25.....	<b>25</b>	56.22.25P.....	<b>42</b>
11.40PIN.....	<b>26</b>	23.76.25PIN.....	<b>27</b>	40.01-2.....	<b>50</b>	56.24.23P.....	<b>44</b>
11.41.23.....	<b>28, 54</b>	23101SYN.....	<b>34</b>	40.01-6.....	<b>50</b>	56.24.25P.....	<b>44</b>
12.01.....	<b>33</b>	23102SYN.....	<b>34</b>	40.02-2.....	<b>50</b>	56.30P.....	<b>44</b>
12.02.....	<b>34</b>	23103SYN.....	<b>34</b>	40.02-6.....	<b>50</b>	56.30.23P.....	<b>44</b>
12.03.....	<b>34</b>	26.01.....	<b>31</b>	41.00.....	<b>53</b>	56.30.25P.....	<b>44</b>
12.03.25.....	<b>34</b>	27.01.....	<b>31</b>	51.20.....	<b>29</b>	56.34.25P.....	<b>44</b>
12.06.25.....	<b>50</b>	28.04.....	<b>30</b>	51.25.....	<b>29</b>	56.44P.....	<b>29, 55</b>
12.06.27.....	<b>50</b>	28.05.....	<b>31</b>	51.35.....	<b>29</b>	56.45.23P.....	<b>29, 55</b>
12.12.....	<b>34</b>	28.06.....	<b>31</b>	51.45.....	<b>29</b>	56.46.23P.....	<b>29, 55</b>
12.13.....	<b>34</b>	28.07.....	<b>31</b>	52.00.....	<b>49</b>	56.50.25P.....	<b>43</b>
12.21.....	<b>26</b>	28.21.....	<b>30</b>	52.20.....	<b>49</b>	56.50.29P.....	<b>43</b>
12.21PIN.....	<b>26</b>	28.25.....	<b>34</b>	53.00.23.....	<b>49, 51</b>	56.54.25P.....	<b>43</b>
12.40.....	<b>35</b>	28.30.....	<b>31, 54</b>	53.00.25.....	<b>49, 51</b>	56.54.27P.....	<b>43</b>
12.41.....	<b>35</b>	28.30.23.....	<b>54</b>	53.00.27.....	<b>49, 51</b>	56.65.....	<b>21</b>
13.01.....	<b>33</b>	29.22.....	<b>28</b>	55.21.....	<b>38</b>	56.65.23.....	<b>21</b>
14.01.....	<b>33</b>	29.22.23.....	<b>28</b>	55.21.23.....	<b>38</b>	56.65.25.....	<b>21</b>
		30.01.....	<b>32</b>	55.21.25S.....	<b>38</b>	56.65.27.....	<b>21</b>
		30.02.....	<b>32</b>	55.21.27.....	<b>38</b>	57.27.....	<b>42</b>
		30.06.....	<b>32</b>	55.23.....	<b>38</b>	57.30.....	<b>42</b>
		30.08.....	<b>32</b>	55.23.23.....	<b>38</b>	57.40.....	<b>42</b>
		30.08.23.....	<b>34</b>	55.23.25.....	<b>38</b>	58.2023-3MM.....	<b>54</b>
		30.08.25.....	<b>34</b>	55.23.27.....	<b>38</b>	58.20-6.....	<b>21</b>
		30.08.27.....	<b>34</b>	55.25.23.....	<b>38</b>	58.23-3MM.....	<b>54</b>
		30.09.....	<b>32</b>	55.25.25.....	<b>38</b>	58.23-3MMV.....	<b>54</b>
		30.10.....	<b>32</b>	55.26.....	<b>38</b>	58.23STHV.....	<b>21</b>
		30.10.23.....	<b>32</b>	55.26E.....	<b>38</b>	58.25STHV.....	<b>21</b>
		30.18.23.....	<b>34</b>	55.26.23.....	<b>38</b>	70020S.....	<b>19</b>



70020SP.....	<b>19</b>	BL5293LC.....	<b>39</b>	E7931.....	<b>51</b>
70023E.....	<b>19</b>	BL5295.....	<b>39</b>	ES5000S.....	<b>18</b>
70023S.....	<b>19</b>	BL5295C.....	<b>39</b>	LP7519.....	<b>17</b>
70023SP.....	<b>19</b>	BL5295L.....	<b>39</b>	LP7520.....	<b>17</b>
70025E.....	<b>19</b>	BL5295LC.....	<b>39</b>	LP7523.....	<b>17</b>
70025S.....	<b>19</b>	BL5320WX.....	<b>15</b>	LP7525.....	<b>17</b>
70025SP.....	<b>19</b>	BL5323WVX.....	<b>15</b>	MVS1038A.....	<b>48</b>
70123B.....	<b>19</b>	BL5325WVX.....	<b>15</b>	MVS1069.....	<b>48</b>
70125B.....	<b>19</b>	BL5420WX.....	<b>12</b>	MVS1070.....	<b>48</b>
70933.....	<b>20</b>	BL5423WVX.....	<b>12</b>	MVS1600A.....	<b>7</b>
70935.....	<b>20</b>	BL5425WVX.....	<b>12</b>	P2-PLAMP.....	<b>42</b>
BL14304.....	<b>16</b>	BL5427MVXB.....	<b>10</b>	PLAMP.....	<b>42</b>
BL14334.....	<b>16</b>	BL5520WX.....	<b>14</b>	S2050B.....	<b>51</b>
BL14455.....	<b>6</b>	BL5523WVX.....	<b>14</b>	S2070C.....	<b>51</b>
BL15455.....	<b>6</b>	BL5525WVX.....	<b>14</b>	SE5420M.....	<b>8</b>
BL2272.....	<b>7</b>	BL5527MVXB.....	<b>11</b>	SE5420W.....	<b>8</b>
BL2295.....	<b>6</b>	BL5620.....	<b>17</b>	SE5423MV.....	<b>8</b>
BL2296.....	<b>6</b>	BL5623S.....	<b>17</b>	SE5423WV.....	<b>8</b>
BL2298.....	<b>7</b>	BL5625S.....	<b>17</b>	SE5425MV.....	<b>8</b>
BL3170.....	<b>7</b>	BL5626.....	<b>17</b>	SE5425MVB.....	<b>10</b>
BL3234.....	<b>7</b>	BL5627S.....	<b>17</b>	SE5425WV.....	<b>8</b>
BL3235.....	<b>7</b>	BL5720.....	<b>45</b>	SE5425WVB.....	<b>10</b>
BL3236.....	<b>7</b>	BL5723.....	<b>45</b>	SE5427MVB.....	<b>10</b>
BL3237.....	<b>7</b>	BL5725.....	<b>45</b>	SE5520M.....	<b>9</b>
BL3239.....	<b>7</b>	BL5820.....	<b>45</b>	SE5520W.....	<b>9</b>
BL3240.....	<b>7</b>	BL5820D.....	<b>45</b>	SE5523MV.....	<b>9</b>
BL3241.....	<b>7</b>	BL5823.....	<b>45</b>	SE5523WV.....	<b>9</b>
BL3242.....	<b>7</b>	BL5823D.....	<b>45</b>	SE5525MV.....	<b>9</b>
BL3244.....	<b>7</b>	BL5823M.....	<b>45</b>	SE5525MVB.....	<b>11</b>
BL3245.....	<b>7</b>	BL5825.....	<b>45</b>	SE5525WV.....	<b>9</b>
BL3246.....	<b>7</b>	BL5825D.....	<b>45</b>	SE5525WVB.....	<b>11</b>
BL3247.....	<b>7</b>	BL5825M.....	<b>45</b>	SE5527MVB.....	<b>11</b>
BL3270.....	<b>7</b>	BL5827M.....	<b>45</b>	SE5623.....	<b>17</b>
BL3320.....	<b>7</b>	BL7610.....	<b>18</b>	SE5625.....	<b>17</b>
BL3321.....	<b>7</b>	BL7620.....	<b>18</b>	SE5625B.....	<b>17</b>
BL3394.....	<b>17</b>	BL7623.....	<b>18</b>	SE5627B.....	<b>17</b>
BL4363.....	<b>6</b>	BL7625.....	<b>18</b>	SUK01.....	<b>31</b>
BL4364.....	<b>7</b>	CX5710.....	<b>18</b>	SUK02.....	<b>31</b>
BL4390.....	<b>6</b>	CX5720.....	<b>18</b>	TR906061.....	<b>18</b>
BL4391.....	<b>6</b>	CX9400.....	<b>51</b>	TR9066.....	<b>18</b>
BL4393.....	<b>6</b>	CX9404.....	<b>51</b>	TR9085.....	<b>18</b>
BL4394.....	<b>6</b>	CX9430.....	<b>51</b>	TR9088.....	<b>18</b>
BL5220WX.....	<b>13</b>	D00.01SL.....	<b>24</b>	TR9089.....	<b>18</b>
BL5223WVX.....	<b>13</b>	D00.01SS.....	<b>24</b>	V30PD.....	<b>53</b>
BL5225WVX.....	<b>13</b>	D09.04.....	<b>25</b>	VBCD.....	<b>53</b>
BL5234.....	<b>45</b>	D09.40.....	<b>25</b>	VFD.....	<b>53</b>
BL5235.....	<b>45</b>	D09.40.25.....	<b>25</b>	VFIHP.....	<b>20</b>
BL5235S.....	<b>45</b>	D23.76.....	<b>27</b>	VFIMOD.....	<b>20</b>
BL5237S.....	<b>45</b>	D23.76.23.....	<b>27</b>	VFLATSSVD.....	<b>53</b>
BL5240.....	<b>16</b>	D23.76.25.....	<b>27</b>	VMD.....	<b>53</b>
BL5241.....	<b>16</b>	D28.30.23.....	<b>54</b>	VVITAHOSE-E.....	<b>19</b>
BL5275.....	<b>7</b>	D29.02.....	<b>28</b>	VVITAHOSE-Y.....	<b>19</b>
BL5277.....	<b>7</b>	D31.01.....	<b>33</b>	VVITFUSE-1.0.....	<b>19</b>
BL5280.....	<b>21</b>	D31.06.....	<b>33</b>	VVITFUSE-2.0.....	<b>19</b>
BL5281.....	<b>16</b>	D31.08.....	<b>33</b>	VVITOME.....	<b>19</b>
BL5282.....	<b>16</b>	D31.08.23.....	<b>33</b>	VWFD.....	<b>53</b>
BL5283.....	<b>16</b>	D31.08.25.....	<b>33</b>		
BL5284.....	<b>16</b>	D39.03.23.....	<b>25</b>		
BL5285.....	<b>16</b>	D39.03.25.....	<b>25</b>		
BL5286.....	<b>16</b>	D39.03.27.....	<b>25</b>		
BL5290.....	<b>39</b>	D39.08.....	<b>25</b>		
BL5290ASP.....	<b>39</b>	D39.08.23.....	<b>25</b>		
BL5290ASPH.....	<b>39</b>	D39.08.25.....	<b>25</b>		
BL5290C.....	<b>39</b>	D39.10.23.....	<b>25</b>		
BL5290L.....	<b>39</b>	D39.10.25.....	<b>25</b>		
BL5290LC.....	<b>39</b>	D39.15.23.....	<b>26</b>		
BL5293.....	<b>39</b>	D39.15.25.....	<b>26</b>		
BL5293ASP.....	<b>39</b>	D4600A.....	<b>21</b>		
BL5293ASPH.....	<b>39</b>	E7924.....	<b>51</b>		
BL5293C.....	<b>39</b>	E7925.....	<b>51</b>		
BL5293L.....	<b>39</b>	E7930.....	<b>51</b>		

# Index

## Vitrectomy

Stellaris Elite™	6
BL15455	6
BL14455	6
Hardware + Accessories	6
BL4363	6
BL2295	6
BL2296	6
BL4390	6
BL4391	6
BL4393	6
BL4394	6
BL3270	7
BL4364	7
BL3170	7
BL3236	7
BL3240	7
BL3245	7
BL3241	7
BL3247	7
BL3237	7
BL3239	7
BL3242	7
BL3320	7
BL3321	7
BL5275	7
BL5277	7
BL3235	7
BL3244	7
BL2272	7
BL2298	7
BL3246	7
BL3234	7
MVS1600A	7
Surgical Packs	8
SE5420M	8
SE5420W	8
SE5423MV	8
SE5423WV	8
SE5425MV	8
SE5425WV	8
SE5520M	9
SE5520W	9
SE5523MV	9
SE5523WV	9
SE5525MV	9
SE5525WV	9
BL5427MVXB	10
SE5425MVB	10
SE5425WVB	10
SE5427MVB	10
BL5527MVXB	11
SE5525MVB	11
SE5525WVB	11
SE5527MVB	11
BL5420WX	12
BL5423WVX	12
BL5425WVX	12
BL5220WX	13
BL5223WVX	13
BL5225WVX	13
BL5520WX	14
BL5523WVX	14

BL5525WVX	14
BL5320WX	15
BL5323WVX	15
BL5325WVX	15
Stellaris®PC	16
BL14304	16
BL14334	16
Valved Trocars	16
BL5281	16
BL5282	16
BL5240	16
BL5283	16
BL5285	16
BL5241	16
BL5284	16
BL5286	16
Vitrectomy Cutters	17
BL5620	17
BL5623S	17
BL5625S	17
BL5626	17
SE5623	17
SE5625	17
BL3394	17
SE5625B	17
SE5627B	17
BL5627S	17
MVR Knives	17
LP7519	17
LP7520	17
LP7523	17
LP7525	17
Vitreoretinal Tamponades	18
TR906061	18
TR9066	18
TR9085	18
TR9088	18
TR9089	18
BL7620	18
BL7623	18
BL7625	18
ES5000S	18
CX5710	18
CX5720	18
BL7610	18
30.90	18
VersaVIT™ 2.0 Packs	19
70020S	19
70023S	19
70025S	19
70020SP	19
70023SP	19
70025SP	19
70023E	19
70025E	19
70123B	19
70125B	19
VersaVIT™ 2.0 Accessories	19
VVITOME	19
VVITAHOSE-Y	19
VVITAHOSE-E	19
VVITFUSE-1.0	19
VVITFUSE-2.0	19
VersaVIT™ Accessories	19
17565-10	20
17570	20
17680	20
17580	20
VFIHP	20
VFIMOD	20
19020	20

19030	20
70933	20
70935	20
30.90	21
56.65	21
56.65.23	21
56.65.25	21
56.65.27	21
58.20-6	21
58.23STHV	21
58.25STHV	21
BL5280	21
D4600A	21

## Instrumentation

Handle Designs	24
D00.01SL	24
D00.01SS	24
ILM	25
39.08PIN	25
39.08.23PIN	25
39.08.25PIN	25
39.08.27PIN	25
D39.08	25
D39.08.23	25
D39.08.25	25
39.28	25
39.28.23	25
39.28.25	25
09.12PIN	25
09.12.23PIN	25
09.12.25PIN	25
09.10.23PIN	25
09.10.25PIN	25
D39.10.23	25
D39.10.25	25
09.10	25
09.10.23	25
09.10.25	25
09.26	25
39.12.23PIN	25
39.12.25PIN	25
D09.40	25
D09.40.25	25
09.50	25
09.50.25	25
11.28	25
11.28.23	25
11.23	25
11.23.23	25
09.04PIN	25
39.03.23PIN	25
39.03.25PIN	25
39.03.27PIN	25
D09.04	25
D39.03.23	25
D39.03.25	25
D39.03.27	25
39.23	25
39.23.23	25
39.23.25	25
39.18.23PIN	26
39.24	26
Heavy Membrane	26
11.04PIN	26
11.04.23PIN	26
11.04.25PIN	26
11.24	26
11.06PIN	26

11.06.23PIN	26
11.06.25PIN	26
11.26	26
Dual Purpose	26
39.07.23PIN	26
39.07.25PIN	26
39.27.23	26
39.27.25	26
39.15.23PIN	26
39.15.25PIN	26
39.15.27PIN	26
D39.15.23	26
D39.15.25	26
Foreign Body	26
11.40PIN	26
11.40	26
11.20PIN	26
11.20	26
12.21PIN	26
12.21	26
Scissors	27
23.76PIN	27
23.76.23PIN	27
23.76.25PIN	27
D23.76	27
D23.76.23	27
D23.76.25	27
23.26	27
23.26.23	27
23.26.25	27
23.27	27
23.56	27
23.56.23	27
23.56.25	27
10.24PIN	27
10.24.23PIN	27
10.24.25PIN	27
10.24	27
10.24.23	27
10.24.25	27
23.29	27
23.29.23	27
10.21	27
Pick Forceps	28
D29.02	28
29.22	28
29.22.23	28
09.05PIN	28
09.05.23PIN	28
09.05.25PIN	28
09.25	28
09.25.23	28
09.25.25	28
11.41.23	28
Subretinal	28
39.21	28
23.23	28
21.21	28
Myopic	28
09.10.23PIN-M	28
09.10.25PIN-M	28
39.15.23PIN-M	28
39.15.25PIN-M	28
Illuminated Disposables	29
56.44P	29
56.45.23P	29
56.46.23P	29
Scleral Markers	29
51.20	29
51.35	29
51.45	29

51.25.....	29
Picks + Spatulas.....	<b>30</b>
20.02.....	30
20.02.23.....	30
20.02.25.....	30
20.41.....	30
20.41.23.....	30
20.41.25.....	30
28.21.....	30
20.31.....	30
28.04.....	30
20.01.....	31
28.30.....	31
20.30.....	31
26.01.....	31
27.01.....	31
28.07.....	31
28.05.....	31
28.06.....	31
SUK01.....	31
SUK02.....	31
Cannulas + Needles.....	<b>32</b>
31.00.....	32
30.10.....	32
30.10.23.....	32
30.01.....	32
30.02.....	32
30.06.....	32
30.08.....	32
30.09.....	32
31.01.....	32
31.02.....	32
02.03.....	32
08.01.....	33
12.01.....	33
13.01.....	33
14.01.....	33
15.01.....	33
D31.01.....	33
D31.06.....	33
D31.08.....	33
D31.08.23.....	33
D31.08.25.....	33
12.02.....	34
12.03.....	34
12.12.....	34
12.13.....	34
12.03.25.....	34
30.08.23.....	34
30.08.25.....	34
30.08.27.....	34
30.18.23.....	34
30.18.25.....	34
30.76.....	34
30.80.25.....	34
19020.....	34
19030.....	34
Instrument Trays.....	<b>34</b>
28.25.....	34
23102SYN.....	34
23101SYN.....	34
23103SYN.....	34
Tip Protectors.....	<b>35</b>
34.71.....	35
34.72.....	35
34.74.....	35
34.78.....	35
Luers + Plugs.....	<b>35</b>
34.80.....	35
12.41.....	35
12.40.....	35

## Lasers Probes

Directional™.....	<b>38</b>
55.26.....	38
55.26.23.....	38
55.26.25.....	38
55.26.25S.....	38
55.26.27.....	38
55.26E.....	38
55.26.23E.....	38
55.26.25E.....	38
55.36E.....	38
55.36.23E.....	38
55.36.25E.....	38
55.27.....	38
55.27.23.....	38
55.27.25.....	38
55.27.27.....	38
Straight.....	<b>38</b>
55.21.....	38
55.21.23.....	38
55.21.25S.....	38
55.21.27.....	38
Curved.....	<b>38</b>
55.25.23.....	38
55.25.25.....	38
55.23.....	38
55.23.23.....	38
55.23.25.....	38
55.23.27.....	38
Illuminated.....	<b>39</b>
55.69P.....	39
55.69.23P.....	39
55.69.25P.....	39
55.62.23P.....	39
55.62.25P.....	39
55.70.23P.....	39
55.70.25P.....	39
55.47P.....	39
55.47.23P.....	39
55.47.25P.....	39
55.48.23P.....	39
55.48.25P.....	39
Stellaris® Disposables.....	<b>39</b>
BL5290.....	39
BL5293.....	39
BL5295.....	39
BL5290C.....	39
BL5293C.....	39
BL5295C.....	39
BL5290L.....	39
BL5293L.....	39
BL5295L.....	39
BL5290LC.....	39
BL5293LC.....	39
BL5295LC.....	39
BL5290ASPH.....	39
BL5293ASPH.....	39
BL5290ASP.....	39
BL5293ASP.....	39

## Illumination

Light Source Accessories.....	<b>42</b>
PLAMP.....	42
P2-PLAMP.....	42
Illumination Adapters.....	<b>42</b>
57.27.....	42
57.30.....	42
57.40.....	42

Endoilluminators.....	<b>42</b>
56.21P.....	42
56.21.23P.....	42
56.21.25PS.....	42
56.12P.....	42
56.12.23P.....	42
56.12.25P.....	42
56.22P.....	42
56.22.23P.....	42
56.22.25P.....	42
56.20P.....	42
56.20.23P.....	42
56.20.25P.....	42
56.02P.....	42
56.02.23P.....	42
56.02.25PS.....	42
56.02.27P.....	42
Chandeliers.....	<b>43</b>
56.54.25P.....	43
56.54.27P.....	43
56.50.25P.....	43
56.50.29P.....	43
Infusion Chandeliers.....	<b>44</b>
56.30P.....	44
56.30.23P.....	44
56.30.25P.....	44
56.34.25P.....	44
Multifunctional.....	<b>44</b>
56.07.23P.....	44
56.07.25PS.....	44
56.14.....	44
56.00.25P.....	44
56.24.23P.....	44
56.24.25P.....	44
Stellaris® Disposables.....	<b>45</b>
BL5270.....	45
BL5273.....	45
BL5275.....	45
BL5823M.....	45
BL5825M.....	45
BL5827M.....	45
BL5820D.....	45
BL5823D.....	45
BL5825D.....	45
BL5820.....	45
BL5823.....	45
BL5825.....	45
BL5235S.....	45
BL5237S.....	45
BL5234.....	45
BL5235.....	45

## Other Products

DDMSTM.....	<b>48</b>
20.04.....	48
20.04.23.....	48
20.04.25S.....	48
20.04.27.....	48
20.07.....	48
20.16.....	48
20.16.23.....	48
20.16.25.....	48
20.18.....	48
Infusion Cannulas.....	<b>48</b>
30.32.23.....	48
30.32.25.....	48
MVS1038A.....	48
MVS1069.....	48
MVS1070.....	48

30.20.....	49
30.21.....	49
30.22.....	49
Bipolar.....	<b>49</b>
53.00.23.....	49
53.00.25.....	49
53.00.27.....	49
Knives.....	<b>49</b>
52.00.....	49
52.20.....	49
Iris Retractors.....	<b>50</b>
40.01-2.....	50
40.01-6.....	50
40.02-2.....	50
40.02-6.....	50
Visco Dissector.....	<b>50</b>
12.06.25.....	50
12.06.27.....	50
Stellaris® Bipolar.....	<b>51</b>
E7930.....	51
E7931.....	51
E7924.....	51
E7925.....	51
CX9404.....	51
CX9400.....	51
CX9430.....	51
S2050B.....	51
S2070C.....	51
53.00.23.....	51
53.00.25.....	51
53.00.27.....	51
I-PACK™ Components.....	<b>52</b>
18060.....	52
18061.....	52
18064.....	52
18067.....	52
18069.....	52
18069W.....	52
Lenses.....	<b>53</b>
41.00.....	53
VFD.....	53
VBCD.....	53
VMD.....	53
VWFD.....	53
V30PD.....	53
VFLATSSVD.....	53
Pediatric.....	<b>54</b>
58.2023-3MM.....	54
58.23-3MM.....	54
58.23-3MMV.....	54
11.41.23.....	54
39.18.23PIN.....	54
28.30.....	54
28.30.23.....	54
D28.30.23.....	54
10.24PIN.....	54
10.24.23PIN.....	54
10.24.25PIN.....	54
56.45.23P.....	55
56.46.23P.....	55
56.14.....	55
56.44P.....	55

800.338.2020  
bauschretina.com  
customerservice@bauschretina.com

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.  
Any other product/brand names and/or logos are trademarks of the respective owners.  
© Bausch & Lomb Incorporated.

SRG.0190.USA.18



WARNING: Class 4 Laser Product—visible laser radiation.  
Avoid eye or skin exposure to direct or scattered radiation.

