

Densiron® 68

"Heavier-than-water" endotamponade with exceptional features

EP 0859751, US 6,626,126 B1

Proven
Quality!



Increased density of 1.06 g/cm³ to efficiently tamponade the inferior retina

Purely physical mix of two intensively researched and CE-certified liquid implants (F6H8® and Siluron® 5000)

Mid-range viscosity of 1400 mPas

Easy application and removal of the endotamponade

High storage stability even at low temperatures (no phase separation up to 2 °C)

Reliable and successful in clinical applications since more than 10 years

G-80910 Densiron® 68 vial

10 ml vial, sterile

G-80920 Densiron® 68 syringe

10 ml syringe, sterile

Geuder®
Precision made in Germany

www.geuder.com

FLUORON®
Leading in purity and variety

Regarding product availability in the individual countries please contact Geuder AG or your distributor!

69045

05GLG5